




Feasibility/Collaboration Proposal: Model calibration as a Tool for Material Characterisation

Bilal Ahmed, Faraj A Shmam,
Neeru Bala, [Daniel Markl](#) , Rachel Smith

University of Sheffield & University of Strathclyde




Agile & Flexible Development & Manufacturing




Agile & Flexible Development & Manufacturing

to reduce speed & cost and increase efficiency of how medicines are developed and manufactured to meet changing demands.




Product Performance

Tablet Disintegration – Dissolution Relationship




Product Performance


Disintegration and Dissolution




Disintegration of Granulated Tablets




Experimental Investigation




Investigating swelling of single particles



Investigating liquid uptake of granules




Investigating liquid uptake and swelling of granules




Single Particle Swelling Characterisation

In-situ monitoring and modelling of single particle swelling




Extraction of characteristic swelling properties and swelling-limiting mechanisms









Single Particle Swelling Characterisation

1 *In-situ* images of single particle swelling




2 Single particle swelling modelling and parameter estimation




PH101/102: microcrystalline cellulose grades
 CCS: croscarmellose sodium
 SSG: sodium starch glycolate
 L-HPC: hydroxypropyl cellulose

D : Diffusion coefficient
 Q^{max} : maximum absorption ratio
 Δr^{max} : swelling capacity


Experimental Investigation




Investigating swelling of single particles




Investigating liquid uptake of granules



Investigating liquid uptake and swelling of granules





1 Granule generation via high shear wet granulator



Material	Function	%
Aspirin	API	10
Microcrystalline cellulose	Filler, partial lubricant & disintegrant	80-90
Sodium Starch Glycolate (SSG)	Superdisintegrant	2-7%
Water/PEG	Liquid binder	-

2 Granule swelling performance assessment




Material	Size, x (mm)	Porosity, ε (%)	L/S	Flow rate, \dot{Q} (mL/min)
MCC	1.9-2.5	47	1	10
MCC+SSG	1.8-1.9	25	1	10
Aspirin+MCC	1.9-2.5	30-62	0.8-1.2	10
Aspirin+MCC +SSG	1.9-2.5	35-52	1	10

Granule Swelling Characterisation


Microcrystalline Cellulose

MCC



- ▶ Surface erosion and slow crack propagation observed over time
- ▶ No excessive granule movement, challenges with bubble formation

MCC + SSG




- ▶ Adding superdisintegrant caused sudden expansion and swelling over time
- ▶ Slow loosening of microstructure and strain recovery observed

Granule Swelling Characterisation


Aspirin and Excipients

Aspirin + MCC




- ▶ Aspirin added to MCC (no superdisintegrant)
- ▶ Erosion and crack propagation observed with slow swelling

Aspirin + MCC + SSG




- ▶ All components included within granule
- ▶ Controlled crack propagation and erosion occurs. Gradual swelling


Experimental Investigation




Investigating swelling of single particles



Investigating liquid uptake of granules




Investigating liquid uptake and swelling of granules





Granule Liquid Absorption and Swelling Characterisation

1 Granule generation via high shear wet granulator



Material	Function	%
Aspirin	API	10
Microcrystalline cellulose	Filler, partial lubricant & disintegrant	80-90
Sodium Starch Glycolate	Superdisintegrant	2-7%
Water/PEG	Liquid binder	-


2 Granule swelling performance assessment



Material	Size, x (mm)	Porosity, ε (%)	L/S	Flow rate, \dot{Q} (mL/min)
MCC	1.9-2.5	-	1	10
MCC+SSG	1.8-1.9	-	1	10
Aspirin+MCC	1.9-2.5	30-62	0.8-1.2	10
Aspirin+MCC +SSG	1.9-2.5	-	1	10

Optical Coherence Tomography for Material Characterisation






- Optical coherence tomography is a non-destructive, contactless and fast method to resolve sub-surface structures.
- It generates 1D, 2D and 3D depth profiles within seconds.
- Acquisition rate of 1D profiles can be as high as 230 kHz (\rightarrow 230 cross-sectional images per second)







Selected Application of OCT

Tablet Coating Thickness Characterisation




Tablet Coating Characterisation



Pellet Coating Characterisation




Disintegration Analysis using Optical Coherence Tomography



OCT conditions:

- Lens: 18 mm focal length
- Burst of images collected upon liquid addition and granule events
- Only 2D images analysed







Product Structure




$$D_p \text{ (mm)} = 2.3$$

$$\varepsilon \text{ (%) } = 47 \quad \text{● MCC}$$

Surface Imaging







Product Structure
Surface Imaging

 D_p (mm) = 2.5

 ϵ (%) = 35

 MCC




 Aspirin

 SSG (2%)

 D_p (mm) = 2.5


 ϵ (%) = 52

 MCC

 Aspirin

 SSG (7%)


Advanced performance tests



1 Identify performance-controlling mechanisms
2 Develop data-driven, mechanistic and/or semi-mechanistic models

